



**OHIO DEPARTMENT OF HEALTH**  
**BUREAU OF ALCOHOL AND DRUG TESTING**

APPROVAL OF INSTRUMENT CHECK SOLUTION, OAC 3701-53-04

**Lot or Batch # 09210 Target Value 0.100 g/210L**

This instrument check solution is to be used in a simulator for the instrument check of approved evidential breath testing instruments, as defined in OAC 3701-53-02.

This instrument check solution contains 1.217 mg/mL ethyl alcohol in purified water. When used according to instrument check forms, it will produce a result at or within  $\pm 0.005$  g/210L of the target value in the approved breath testing instrument, when that instrument is in proper working condition.

This instrument check solution shall not be used more than three months after its date of first use or more than one year after its date of manufacture, July 20, 2009. Whichever comes first.

**This solution expires July 20, 2010**

Approved this 13<sup>th</sup> day of August, 2009

A handwritten signature in black ink, appearing to read "Timothy J. Fackler".

Director of Health