



Vaccine Storage Guidelines



COMBINATION REFRIGERATOR/FREEZER

- ✓ Store vaccine in original packaging; group by type and clearly label the designated space for each vaccine.
- ✓ Keep vaccine 2-3 inches away from walls and other boxes.
- ✓ Post "Do Not Unplug" stickers on outlets. Plug in ONLY 1 unit per outlet.
- ✓ Place thermometer probe in the center of the unit. Post a temperature log on the door.
- ✓ Never store food in the refrigerator or freezer.



For information about vaccine storage and handling call the Immunization Program:
1-800-282-0546

----- Dashed lines show usable space.

